

STIMULATION WITH PBMC PULSED WITH KNOWN SPECIFIC CTL EPITOPES Method 2 (of 2)

1. PURPOSE AND DESCRIPTION

- 1.1 These methods measure precursor frequencies recognizing a specific HIV-1 gene product (method 1) or CTL epitope (method 2) by limiting dilution assay.
- 1.2 The data are useful in determining quantitative changes in HIVspecific effector responses in a given individual over time. This should provide a better estimate of the effects of a given intervention on CTL responses than qualitative measurements of specific lysis to a given gene product.

2. SPECIMEN

- 2.1 Anticoagulated blood (ACD or heparin) fresh blood, processed within 24 hours of collection, can be shipped overnight at room temperature (Preliminary)

3. MATERIALS AND REAGENTS

- 3.1 20 mL whole heparinized blood
- 3.2 RPMI (Bio-Whittaker, Walkersville, MD)
- 3.3 RPMI + 10% FBS (Biocell, Rancho Dominguez, CA; lot #2262) + penicillin 100 U/mL; streptomycin 100 mg/mL and L-glutamine 1 mM/mL (R10)
- 3.4 recombinant vaccinia viruses: vLac (vSC-8) and vDK1 (HIV-1 gag, D. Kouritzkes)
- 3.5 HIV peptides: SLYN, ILKE, KIRL, QVPL, TLYC, AIPQ, TPGP, GGKK, GPKV, PLKE, NSSK, PPIP, and VPLR (QCB, Inc., Hopkinton, MA)
- 3.6 Lymphocult T (Biotest, Dreieich, Germany)
- 3.7 96-well round-bottom microtiter plates (Costar, Cambridge, MA)

4. INSTRUMENTATION

- 4.1 Radioactive isotope counter capable of counting Cr release, calibration 51 according to manufacturer's directions.

5. PROCEDURE

- 5.1 Day 0
 - 5.1.1 Isolate PBMC from whole blood using Ficoll-Hypaque density centrifugation:
 - 5.1.1.1 Layer whole blood on Ficoll; spin at 2000 rpm for 20 min. with no

brake.

- 5.1.1.2 Wash once with PBS + 4 mM EDTA at 1500 rpm for 10 min with brake
- 5.1.1.3 Wash second time with PBS + 4 mM EDTA at 1200 rpm for 10 min with brake.

5.1.2 Preparation of Stimulator Cells

- 5.1.2.1 Take 2×10^6 cells from PBMCs isolated from whole blood for 6 stimulator cells; wash by centrifugation at 1500 rpm for 5 minutes.
- 5.1.2.2 Discard supernatant, resuspend cell pellet in 500 ml working peptide solution (100 ml peptide/mL)
- 5.1.2.3 Incubate for 2 h at 37°C; 5% CO₂
- 5.1.2.4 Irradiate cells at 3000 rad
- 5.1.2.5 Add 2 mL 100 U/mL lymphocult T + 8 mL R10
- 5.1.2.6 Add 50 ml/well in 96-well plates (50,000 cells/well)

5.1.3 Preparation of Effector Cells

- 5.1.3.1 Prepare 3 mL of 1.6×10^6 cells/mL in R10 5
- 5.1.3.2 Add 50 ml R10 to rows "4000 cells/well" through "63 cells/well"
- 5.1.3.3 Add 50 mL of cell suspension to rows "8000 cells/well" and "4000 cells/well"
- 5.1.3.4 Take 50 mL from "4000 cells/well" and add to "2000 cells/well." Repeat through "63 cells/well." See Figure 1 for plating sample.
- 5.1.3.5 Final volume is 100 ml (including stimulator cells) in all wells; incubate at 37°C; 5% CO₂
- 5.1.3.6 Microcultures are fed at 5 and 10 days by the addition (no removal of media) of 50 ml lymphocult T-containing media (1 mL 100 'U'/mL lymphocult T + 9 mL R10). Total volume after day 10=200 ml.

5.2 CTL assay

- 5.2.1 After 14 days of culture, the effector cells are assessed for cytotoxicity against Cr-labeled target cells. 51
- 5.2.2 The evening before the addition of targets (Day 13) the effector cells are mixed in their individual wells and equal volumes are transferred to 96-well round bottomed plates (depending on the number of targets; each target is tested against the entire plate: i.e. if 4 targets, then each well of effectors is split to 5 new wells in 5 new 96-well plates to include control). See Figure 2 for plate setup.
- 5.2.3 Add 50 ml R10 to each well of target 96-well plate.
- 5.2.4 Add 100 ml of appropriate target cells to each well (105/ml).
- 5.2.5 Add 100 ml of target cells to 24 wells of separate Min/Max plate (24 wells for each target begin tested).
- 5.2.6 Add 100 ml R10 to the first 12 wells for each target in the Min/Max

plate. Add 100 ml 5% Triton solution to the final 12 wells for each target. See Figure 3 for plate setup.

5.2.7 Spin plates at 500 g for 2 min with no brake.

5.2.8 Incubate the plates for 5 hours at 37 C, 5% CO₂.

5.2.9 Spin plates again before collecting the supernatants.

5.2.10 After 5 hours, harvest 30 ml supernatants, and count Cr. 51

5.3 Preparation of Targets

5.3.1 For each target to be tested: centrifuge 2.5×10^6 autologous BLCL in a 15cc tube, discard the supernatant, and resuspend cells in 200 ml R10.

5.3.2 Add 100 mCi 51Cr.

5.3.3 Add 200 ml 100 mg/ml working solution peptide (for peptide-infected targets; 100 ml/ml)

5.3.4 Infect with recombinant vaccinia (MOI=1; for vaccinia-infected targets)

5.3.5 Incubate at 37°C in 5% CO₂ overnight (12 ~ 16h).

5.3.6 Wash target cells by centrifugation 3X

5.3.7 Resuspend cells at 10^5 /ml

5.3.8 Add 100 ml of target cells to each well

- e.g. 1 1' plate control peptide
- 2 2' plate TLYC (p17)
- 3 3' plate AIPQ (RT)
- 4 4' plate vLac (control construct)
- 5 5' plate vDK-1 (gag)

6. CALCULATIONS

6.1 Determine the number of negative wells for each concentration of responder cells and for each target. A positive response is arbitrarily defined as greater than 3 standard deviations above the mean spontaneous release.

6.2 Plot the fraction of negative wells (y axis, log scale) for each cell concentration (x axis). The frequency of antigen-specific pCTL is estimated by the cell concentration in which 37% of the wells are negative for lysis. This is calculated using a chi-square minimization method (Ref: Clouse et al. 1989, J Clin Microbiol 27: 2316-2323). (Program for these calculations can be provided on disk upon request).

7. FREQUENCY/TOLERANCE/CORRECTIONS

Studies preliminary, this information pending

8. EXPECTED VALUES

pCTL frequencies varying from 0-1000 per 106 PBMC

9. PROCEDURE NOTES

Pending, analysis too preliminary

10. LIMITATIONS OF METHOD

- 10.1 Lack of recognition of epitopes included in peptide panel in method 2.
- 10.2 Sensitivity of method 1.
- 10.3 Need for lab workers to be immune to vaccinia (Method 1)
- 10.4 Lysis of control targets excessive (preexisting vaccinia or EBV responsiveness)
- 10.5 Viability and quantity of EBV-transformed B-LCL sufficient to perform the assay.

11. METHOD VALIDATION: INTRA ASSAY VARIABILITY AND INTRA-SUBJECT VARIABILITY

- 11.1 Work in progress for these determinations.

12. REFERENCES

To be provided

13. EFFECTIVE DATE/SCHEDULE

9/30/96

14. AUTHORS

- 14.1 M. Juliana McElrath, M.D., Ph.D.
- 14.2 Norman Letvin, M.D.

Procedure: ACTG Lab Man Cytotoxic T-Lymphocyte Method 2

Prepared by: ACTG Laboratory Technologist Committee

Preparation Date: 01 June 2004

Date Implemented into the Laboratory: _____

Updated on:

Reviewed by:

Date:

Supersedes Archived Manual: DAIDS Virology Manual for HIV Laboratories, Version January 1997